

4-gVirtualXRay_vs_DR-Lungman

January 26, 2023

```
[1]: from IPython.display import display  
from IPython.display import Image
```

```
[2]: from utils import * # Code shared across more than one notebook
```

```
[3]: output_path = "4-output_data/"  
  
if not os.path.exists(output_path):  
    os.mkdir(output_path)
```

1 gVirtualXray vs DR

Main contributors: T. Wen, J. Pointon, J. Tugwell-Allsup and F. P. Vidal

Purpose: We aim to reproduce a real digital radiograph taken with a clinically utilised X-ray equipment.

Material and Methods:

Lungman phantom	3D model of the phantom, source and detector
-----------------	--




1. The CT of a chest phantom has been generated from a real scanner ahead of time.
2. Structures in the reference CT have been segmented and labelled.
3. The resultant surfaces from the segmentations form a virtual lungman model.
4. A digital radiograph was taken with a clinically utilised X-ray equipment.
5. We extract acquisition parameters from the DICOM file and initialise the X-ray simulation

parameters.

6. Using a multi-objective optimisation algorithm, the virtual phantom is registered so that its simulated radiograph closely match the real digital radiograph.

Results: The **zero-mean normalised cross-correlation** is **98.89%**. The **Structural Similarity Index (SSIM)** is **0.94**. The **mean absolute percentage error (MAPE)** is **1.57%**. These results show that the two images are comparable.



The calculations were performed on the following platform:

[4] : `printSystemInfo()`

OS:

Linux 5.3.18-150300.59.54-default
x86_64

CPU:

AMD Ryzen 7 3800XT 8-Core Processor

RAM:

63 GB

GPU:

Name: NVIDIA GeForce RTX 2080 Ti

Drivers: 525.60.13

Video memory: 11 GB

1.1 Import packages

[5] : `%matplotlib inline`

```
import os # Locate files
from time import sleep

import datetime
import math
import numpy as np # Who does not use Numpy?
```

```

import pandas as pd # Load/Write CSV files

# from scipy.stats import pearsonr # Compute the correlation coefficient
# from skimage.util import compare_images # Checkboard comparison between two
# images

# from skimage.util import compare_images # Checkboard comparison between two
# images
from skimage.metrics import structural_similarity as ssim
from sklearn.metrics import mean_absolute_percentage_error as mape
from sklearn.metrics import mean_absolute_error, mean_squared_error
from skimage.metrics import normalized_mutual_information
# from sklearn.metrics.cluster import normalized_mutual_info_score
from skimage.filters import gaussian # Implementing the image sharpening filter
from sklearn.metrics.cluster import normalized_mutual_info_score


import cv2

from tifffile import imread, imwrite # Load/Write TIFF files

# import viewscad # Use OpenSCAD to create STL files

# import pyvista as pv # 3D visualisation
# from pyvista import themes

# import cma # Optimise the parameters of the noise model

# import random
import base64

# import urllib, gzip # To download the phantom data, and extract the
# corresponding Z file

import spekpy as sp # Generate a beam spectrum
from scipy import signal # Resampling the beam spectrum

import SimpleITK as sitk

from gvxrPython3 import gvxr # Simulate X-ray images
gvxr.useLogFile()
from gvxrPython3 import json2gvxr # Set gVirtualXRay and the simulation up
from gvxrPython3.utils import visualise

from scipy.spatial import distance # Euclidean distance

```

```

# Ignore some warnings
import warnings
warnings.filterwarnings("ignore", category=UserWarning)

# Visualisation

import matplotlib
# old_backend = matplotlib.get_backend()
# matplotlib.use("Agg") # Prevent showing stuff

from matplotlib.cm import get_cmap
import matplotlib.pyplot as plt # Plotting
from matplotlib.colors import LogNorm # Look up table
import matplotlib.colors as mcolors

# font = {'family' : 'serif',
#          #'weight' : 'bold',
#          'size'   : 22
#          }
# matplotlib.rc('font', **font)
# matplotlib.rc('text', usetex=True)

import plotly.express as px
# import plotly.graph_objects as go
from plotly.subplots import make_subplots

from pymoo.core.problem import Problem
from pymoo.optimize import minimize
from pymoo.algorithms.soo.nonconvex.cmaes import CMAES

```

SimpleGVXR 2.0.4 (2023-01-18T09:50:01) [Compiler: GNU g++] on Linux
gVirtualXRay core library (gvxr) 2.0.4 (2023-01-18T09:50:00) [Compiler: GNU g++]
on Linux

```
[6]: def standardisation(img):
      return (img - img.mean()) / img.std()
```

2 Preparation of the ground truth image

2.1 Read the real X-ray radiograph from a DICOM file

```
[7]: reader = sitk.ImageFileReader()
reader.SetImageIO("GDCMImageIO")
reader.SetFileName("lungman_data/CD3/DICOM/ST000000/SE000000/DX000000")
reader.LoadPrivateTagsOn()
reader.ReadImageInformation()
volume = reader.Execute()
raw_reference_before_cropping = sitk.GetArrayFromImage(volume)[0]
raw_reference_before_cropping.shape = raw_reference_before_cropping.shape

y_min_id = 200
y_max_id = raw_reference_before_cropping.shape[0] - 170
x_min_id = 50
x_max_id = raw_reference_before_cropping.shape[1] - 100
raw_reference = raw_reference_before_cropping[y_min_id:y_max_id, x_min_id:
    ↪x_max_id]

print("The shape was", raw_reference_before_cropping.shape, "| now it is", ↪
    ↪raw_reference.shape)
```

The shape was (1881, 1871) | now it is (1511, 1721)

```
[8]: imwrite(output_path + '/real_projection-lungman.tif', raw_reference.astype(np.
    ↪single))
```

3 Extract information about the “sizes” of pixels from the DICOM file

The descriptions have been extracted from <https://dicom.innolitics.com>

3.1 Imager Pixel Spacing Attribute

Physical distance measured at the front plane of the Image Receptor housing between the center of each pixel. Specified by a numeric pair - row spacing value (delimiter) column spacing value - in mm.

In the case of CR, the front plane is defined to be the external surface of the CR plate closest to the patient and radiation source.

3.2 Detector Element Physical Size Attribute

Physical dimensions of each detector element that comprises the detector matrix, in mm.

Expressed as row dimension followed by column.

Note:

This may not be the same as Detector Element Spacing (0018,7022) due to the presence of spacing material between detector elements.

3.3 Detector Element Spacing Attribute

Physical distance between the center of each detector element, specified by a numeric pair - row spacing value(delimiter) column spacing value in mm. See Section 10.7.1.3 for further explanation of the value order.

Note:

This may not be the same as the Imager Pixel Spacing (0018,1164), and should not be assumed to describe the stored image.

```
[9]: imager_pixel_spacing = np.array(volume.GetMetaData("0018|1164").split("\\")) .  
    ↪astype(np.single)  
detector_element_physical_size = np.array(volume.GetMetaData("0018|7020")) .  
    ↪split("\\")) .astype(np.single)  
detector_element_spacing = np.array(volume.GetMetaData("0018|7022").split("\\")) .  
    ↪astype(np.single)  
  
print("Imager Pixel Spacing (in mm): ", imager_pixel_spacing)  
print("Detector Element Physical Size (in mm): ", detector_element_physical_size)  
print("Detector Element Spacing (in mm): ", detector_element_spacing)
```


Imager Pixel Spacing (in mm): [0.194556 0.194556]

Detector Element Physical Size (in mm): [0.2 0.2]

Detector Element Spacing (in mm): [0.2 0.2]

We plot the image using a linear look-up table and a power-law normalisation.

```
[10]: displayLinearPowerScales(raw_reference,  
                           "Reference image from the Lungman",  
                           output_path + "/reference-lungman-proj",  
                           log=False,  
                           vmin=-93, vmax=89)
```



Apply a zero-mean, unit-variance normalisation

```
[11]: ground_truth = raw_reference
normalised_ground_truth = standardisation(ground_truth)
imwrite(output_path + '/lungman-normalised_ground_truth.tif', ▾
        normalised_ground_truth.astype(np.single))
```

```
[12]: # meta_data_keys = volume.GetMetaDataKeys()

# print("DICOM fields:")
# for key in meta_data_keys:
#     print(key, volume.GetMetaData(key))
```

```
[13]: # def cleanTags(raw_string):
#     regular_expression = re.compile('<.*?>')
#     clean_text = re.sub(regular_expression, '', raw_string)
#     return clean_text
```

```
[14]: # field = volume.GetMetaData("0033/1022")

# for item in field.split("\n"):
#     if "KV" in item:
#         kv = float(cleanTags(item))

# print(kv)
```

3.4 Extract the image size and pixel spacing using SimpleITK

It will be useful to set the X-ray detector parameters for the simulation, and to display the images in millimetres. We can compare these values with the ones extracted above from the DICOM attributes.

```
[15]: sitk_pixel_spacing = volume.GetSpacing()[0:2]
size = volume.GetSize()[0:2]

print("SITK Pixel Spacing (in mm): ", sitk_pixel_spacing)
print("Image size (in pixels): ", str(size[0]) + " x " + str(size[1]))
```

```
SITK Pixel Spacing (in mm): (0.194556, 0.194556)
Image size (in pixels): 1871 x 1881
```

3.5 Extract the kVp from the DICOM file

It will be useful to generate a realistic beam spectrum.

```
[16]: kVp = float(volume.GetMetaData("0018|0060"))
print("Peak kilo voltage output of the x-ray generator used: ", kVp)
```

```
Peak kilo voltage output of the x-ray generator used: 90.0
```

4 Extract the mAs from the DICOM file

```
[17]: mAs = volume.GetMetaData("0018|1152")
print("The exposure expressed in mAs, for example calculated from Exposure Time and X-ray Tube Current:", mAs)
```

```
The exposure expressed in mAs, for example calculated from Exposure Time and X-ray Tube Current: 1
```

5 Extract the time of X-Ray exposure from the DICOM file

```
[18]: exposure = volume.GetMetaData("0018|1150")
print("Time of X-Ray exposure in msec.:", exposure)
```

```
Time of X-Ray exposure in msec.: 5
```

6 Initialise gVirtualXRay

6.1 Set the experimental parameters (e.g. source and detector positions, etc.)

We use known parameters as much as possible, for example we know the size and composition of the sample. Some parameters are extracted from the DICOM file, such as detector size, pixel resolution, and voltage of the X-ray tube.

```
[19]: distance_source_to_detector = float(volume.GetMetaData("0018|1110"))
distance_source_to_patient = float(volume.GetMetaData("0018|1111"))

print("Distance Source to Detector: ", distance_source_to_detector, "mm")
print("Distance Source to Patient: ", distance_source_to_patient, "mm")

window_size = [800, 450]
# source_position = [0.0, 0.0, source_detector_distance_in_cm - block_thickness_in_cm / 2, "mm"]
# detector_position = [0.0, 0.0, - block_thickness_in_cm / 2, "cm"]
detector_up = [0, 1, 0]
```

Distance Source to Detector: 1800.0 mm
Distance Source to Patient: 1751.0 mm

6.2 Initialise the simulation engine

```
[20]: # Create an OpenGL context
print("Create an OpenGL context:",
      str(window_size[0]) + "x" + str(window_size[1]))
)

gvxr.createWindow(-1, True, "OPENGL")

gvxr.setWindowSize(
    window_size[0],
    window_size[1]
)
```

Create an OpenGL context: 800x450

6.3 Load the scanned object

```
[21]: json2gvxr.initSamples("notebook-4.json", verbose=0)
```

```
[22]: number_of_triangles = 0
```

```

for sample in json2gvxr.params["Samples"]:
    label = sample["Label"]
    number_of_triangles_in_mesh = gvxr.getNumberOfPrimitives(label)
    number_of_triangles += number_of_triangles_in_mesh

```

[23]: skin_bbox = gvxr.getNodeOnlyBoundingBox("Skin", "mm")
print(skin_bbox)

(-159.375, -117.5, -148.39999389648438, 159.375, 107.5, 148.39999389648438)

6.4 Set the source position

[24]: # Set up the beam
print("Set up the beam")
print("\tSource position:", (skin_bbox[0] + skin_bbox[3]) / 2, □
 ↪distance_source_to_detector + skin_bbox[1] + distance_source_to_patient - □
 ↪distance_source_to_detector, (skin_bbox[2] + skin_bbox[5]) / 2, "mm")

gvxr.setSourcePosition((skin_bbox[0] + skin_bbox[3]) / 2, □
 ↪distance_source_to_detector + skin_bbox[1] + distance_source_to_patient - □
 ↪distance_source_to_detector, (skin_bbox[2] + skin_bbox[5]) / 2, "mm")

gvxr.usePointSource()

Set up the beam
Source position: 0.0 1633.5 0.0 mm

6.5 Get the spectrum from the DICOM file

[25]: spectrum = {};
filter_material = "Al" # See email Mon 05/07/2021 15:29
filter_thickness_in_mm = 3 # See email Mon 05/07/2021 15:29

s = sp.Spek(kvp=kVp)
s.filter(filter_material, filter_thickness_in_mm) # Filter by 3 mm of Al
unit = "keV"
k, f = s.get_spectrum(edges=True) # Get the spectrum

min_energy = sys.float_info.max
max_energy = -sys.float_info.max

for energy, count in zip(k, f):
 count = round(count)

 if count > 0:

```

max_energy = max(max_energy, energy)
min_energy = min(min_energy, energy)

if energy in spectrum.keys():
    spectrum[energy] += count
else:
    spectrum[energy] = count

```

Reformat the data

```
[26]: # get the integral nb of photons
nbphotons=0.
energy1 = -1.
energy2 = -1.

for energy in spectrum.keys():

    if energy1<0:
        energy1 = float(energy)
    elif energy2<0:
        energy2 = float(energy)
    nbphotons += float(spectrum[energy])
sampling = (energy2-energy1)

# get spectrum
data = []
for energy in spectrum.keys():
    source = [float(energy),float(spectrum[energy])/(nbphotons*sampling)]
    data.append(source)

data_array = np.array(data)

energies, counts = data_array.T
```


Resample the data to reduce the number of bins

```
[27]: temp_count_set = signal.decimate(counts, 6)
number_of_energy_bins = temp_count_set.shape[0]
temp_energy_set = np.linspace(energies.min(), energies.max(), ↴
    ↪number_of_energy_bins, endpoint=True)

test = temp_count_set > 0
count_set = temp_count_set[test]
energy_set = temp_energy_set[test]
```

Plot the beam spectrum from spekpy and the resampled version

```
[28]: plt.figure(figsize= (20,10))
plt.bar(energy_set, count_set, width=2.8, label="Resampled")
plt.plot(energies, counts, label="Original", color="gray")
plt.xlabel('energy in keV')
plt.ylabel('probability distribution of photons per keV')
plt.legend()
plt.savefig(output_path + "/lungman-projection-spectrum.pdf")
```



6.6 Load the beam spectrum in the simulator

```
[29]: gvxr.resetBeamSpectrum() # To be on the safe side when debugging
for energy, count in zip(energy_set, count_set):
    gvxr.addEnergyBinToSpectrum(energy, unit, count);
```

6.7 Set the X-ray detector

```
[30]: # Set up the detector
print("Set up the detector");
print("\tDetector position:", (skin_bbox[0] + skin_bbox[3]) / 2, skin_bbox[1] +_
    →distance_source_to_patient - distance_source_to_detector, (skin_bbox[2] +_
    →skin_bbox[5]) / 2, "mm")
gvxr.setDetectorPosition((skin_bbox[0] + skin_bbox[3]) / 2, skin_bbox[1] +_
    →distance_source_to_patient - distance_source_to_detector, (skin_bbox[2] +_
    →skin_bbox[5]) / 2, "mm");
```

```
print("\tDetector up vector:", [0, 0, 1])
gvxr.setDetectorUpVector(0, 0, 1);
```

Set up the detector
Detector position: 0.0 -166.5 0.0 mm
Detector up vector: [0, 0, 1]

```
[31]: print("\tDetector number of pixels:", size)
gvxr.setDetectorNumberOfPixels(
    size[0],
    size[1]
);

print("\tPixel spacing:", imager_pixel_spacing)
gvxr.setDetectorPixelSize(
    imager_pixel_spacing[0],
    imager_pixel_spacing[1],
    "mm"
);
```

Detector number of pixels: (1871, 1881)
Pixel spacing: [0.194556 0.194556]

Load the detector response in energy

```
[32]: gvxr.clearDetectorEnergyResponse() # To be on the safe side
gvxr.loadDetectorEnergyResponse("Gate_data/responseDetector.txt",
                                "MeV")
```

6.8 Take a screenshot of the 3D environment

```
[33]: gvxr.displayScene()

gvxr.useNegative()
gvxr.useLighing()
gvxr.useWireframe()
gvxr.setSceneRotationMatrix([0.43813619017601013, 0.09238918125629425, -0.
    ↪8941444158554077, 0.0,
                                0.06627026945352554, 0.9886708855628967, 0.
    ↪13463231921195984, 0.0,
                                0.8964602947235107, -0.11824299395084381, 0.
    ↪4270564019680023, 0.0,
                                0.0, 0.0, 0.0, 1.0])
gvxr.setZoom(2639.6787109375)


gvxr.displayScene()
```

```
[34]: screenshot = gvxr.takeScreenshot()
```

```
[35]: plt.figure(figsize= (10,10))
plt.title("Screenshot")
plt.imshow(screenshot)
plt.axis('off')

plt.tight_layout()

plt.savefig(output_path + '/lungman-projection-screenshot-beam-off.pdf')
plt.savefig(output_path + '/lungman-projection-screenshot-beam-off.png')
```



```
[36]: gvxr.computeXRayImage()
gvxr.displayScene()
```


```
[37]: screenshot = gvxr.takeScreenshot()

plt.figure(figsize= (10,10))
plt.title("Screenshot")
plt.imshow(screenshot)
plt.axis('off')

plt.tight_layout()

plt.savefig(output_path + '/lungman-projection-screenshot-beam-on.pdf')
```

```
plt.savefig(output_path + '/lungman-projection-screenshot-beam-on.png')
```



6.9 Simulation with the default values

```
[38]: # Backup the transformation matrix
global_matrix_backup = gvxr.getSceneTransformationMatrix()
```

```
[39]: def getXRayImage():
    global total_energy_in_MeV

    # Compute the X-ray image
    xray_image = np.array(gvxr.computeXRayImage())

    # Flat-field
    # xray_image /= total_energy_in_MeV

    # Negative
    # x_ray_image = 1.0 - x_ray_image
    return np.flip(xray_image) #np.ones(xray_image.shape).astype(np.single) - xray_image
```

```
[40]: # gvxr.enableArtefactFilteringOnCPU()
gvxr.enableArtefactFilteringOnGPU()
```

```
# gxvr.disableArtefactFiltering() # Spere inserts are missing with GPU
→integration when a outer surface is used for the matrix
```

```
[41]: xray_image = getXRayImage()
```

```
[42]: # total_energy_in_keV = 0.0
# for energy, count in zip(energy_set, count_set):
#     effective_energy = find_nearest(detector_response[:,0], energy / 1000,
→detector_response[:,1])

#     total_energy_in_keV += effective_energy * count

total_energy_in_MeV = gvxr.getTotalEnergyWithDetectorResponse()
```

```
[43]: gvxr.displayScene()
gvxr.useNegative()

gvxr.setZoom(1339.6787109375)
gvxr.setSceneRotationMatrix([0.8227577805519104, 0.1368587613105774, -0.
→5516625642776489, 0.0, -0.5680444240570068, 0.23148967325687408, -0.
→7897683382034302, 0.0, 0.01961756870150566, 0.9631487131118774, 0.
→26820749044418335, 0.0, 0.0, 0.0, 0.0, 1.0])

gvxr.setWindowBackGroundColour(0.5, 0.5, 0.5)

gvxr.displayScene()
```

```
[44]: # gxvr.renderLoop()
```

```
[45]: # print(gvxr.getZoom())
# print(gvxr.setSceneRotationMatrix())
```


```
[46]: screenshot = (255 * np.array(gvxr.takeScreenshot())).astype(np.uint8)
```

```
[47]: plt.figure(figsize= (10,10))
plt.title("Screenshot")
plt.imshow(screenshot)
plt.axis('off')

plt.tight_layout()


plt.savefig(output_path + '/default-screenshot-beam-on-lungman.pdf')
plt.savefig(output_path + '/default-screenshot-beam-on-lungman.png')
```

Screenshot



```
[48]: def logImage(xray_image: np.array, min_val: float, max_val: float) -> np.array:
    log_epsilon = 1.0e-9
    shift_filter = -math.log(min_val + log_epsilon)
    if min_val != max_val:
        scale_filter = 1.0 / (math.log(max_val + log_epsilon) - math.log(min_val + log_epsilon))
    else:
        scale_filter = 1.0
    corrected_image = np.log(xray_image + log_epsilon)
    return (corrected_image + shift_filter) * scale_filter
```

```
[49]: xray_image_cropped = xray_image[y_min_id:y_max_id, x_min_id:x_max_id]
displayLinearPowerScales(1 - logImage(xray_image_cropped, xray_image_cropped.
    .min(), xray_image_cropped.max()),
    "Image simulated using gVirtualXRay before the",
    "registration",
    output_path + "/",
    "gVirtualXRay-before_registration-lungman",
    log=False)
```



6.10 Registration

```
[50]: roi_ground_truth_min = ground_truth.min()
roi_ground_truth_max = ground_truth.max()
standardised_roi_ground_truth = standardisation(ground_truth)

imwrite(output_path + '/standardised_roi_ground_truth-lungman.tif', ↴
        standardised_roi_ground_truth.astype(np.single))
```

Keep a backup of the original source and detector position for registration

```
[51]: source_position_bak = gvxr.getSourcePosition("cm")
detector_position_bak = gvxr.getDetectorPosition("cm")

x_init = [
    # Orientation of the sample
    0.0, 0.0, 0.0,
    # Position of the source
    source_position_bak[0],
    source_position_bak[1],
    source_position_bak[2],
    # Position of the detector
]
```

```

detector_position_bak[0],
detector_position_bak[1],
detector_position_bak[2] #,

# Orientation of the detector
#      det_rotation_angle1 = x[8]
#      det_rotation_angle2 = x[9]

#      1.0 / 3.0, # c1
#      1.0, # gain1
#      0.0, # bias1

#      1.0 / 3.0, # c2
#      1.0, # gain2
#      0.0#, # bias2

#      1.0 / 3.0, # c3
#      1.0, # gain3
#      0.0, # bias3
#      2.0 # gamma
]

pos_offset = 10
angle_offset = 4

x1 = [
    -angle_offset, -angle_offset, -angle_offset,
    source_position_bak[0] - pos_offset, source_position_bak[1] - pos_offset, source_position_bak[2] - pos_offset,
    detector_position_bak[0] - pos_offset, detector_position_bak[1] - pos_offset, detector_position_bak[2] - pos_offset#,
    #      -90, -90,
    #      -10.0,
    #      -10.0,
    #      -10.0,
    #
    #      -10.0,
    #      0.0,
    #      0.0#,
    #
    #      -10.0,
    #      -10.0,
    #      -10.0,
    #      0.0
]

```

```

xu = [
    angle_offset, angle_offset, angle_offset,
    source_position_bak[0] + pos_offset, source_position_bak[1] + pos_offset,
    source_position_bak[2] + pos_offset,
    detector_position_bak[0] + pos_offset, detector_position_bak[1] + pos_offset,
    detector_position_bak[2] + pos_offset #,
    # 90, 90,
    # 10.0,
    # 10.0,
    # 10.0,
    #
    # 10.0,
    # 10.0,
    # 10.0#,
    #
    # 10.0,
    # 10.0,
    # 10.0,
    # 100.0
]

def resetToDefaultParameters():
    gvxr.setSourcePosition(source_position_bak[0], source_position_bak[1],
    source_position_bak[2], "cm")
    gvxr.setDetectorPosition(detector_position_bak[0], detector_position_bak[1],
    detector_position_bak[2], "cm")

    # Restore the transformation matrix
    gvxr.setSceneTransformationMatrix(global_matrix_backup)

```

```

[52]: def setTransformations(x):
    # Orientation of the sample
    sample_rotation_angle1 = x[0]
    sample_rotation_angle2 = x[1]
    sample_rotation_angle3 = x[2]

    gvxr.rotateScene(sample_rotation_angle1, 1, 0, 0)
    gvxr.rotateScene(sample_rotation_angle2, 0, 1, 0)
    gvxr.rotateScene(sample_rotation_angle3, 0, 0, 1)

    # Position of the source
    source_position_x = x[3]
    source_position_y = x[4]
    source_position_z = x[5]

    gvxr.setSourcePosition(

```

```

        source_position_x,
        source_position_y,
        source_position_z,
        "cm"
    )

# Position of the detector
det_position_x = x[6]
det_position_y = x[7]
det_position_z = x[8]

gvxr.setDetectorPosition(
    det_position_x,
    det_position_y,
    det_position_z,
    "cm"
)

# Orientation of the detector
#    det_rotation_angle1 = x[8]
#    det_rotation_angle2 = x[9]

```

```
[53]: def updateXRayImage(x, restore_transformation=True):

    # Backup the transformation matrix
    if restore_transformation:
        matrix_backup = gvxr.getSceneTransformationMatrix()

    # Set the transformations
    setTransformations(x)

    # Compute the X-ray image
    xray_image = getXRayImage()

    #    gvxr.displayScene()
    #    screenshot = gvxr.takeScreenshot()

    # Restore the transformation matrix
    if restore_transformation:
        gvxr.setSceneTransformationMatrix(matrix_backup)

    return xray_image #, screenshot
```

```
[54]: def applyLogScaleAndNegative(image: np.array) -> np.array:
    temp = logImage(image, image.min(), image.max())
    return 1.0 - temp
```

7 Refine the results with CMA-ES

```
[55]: def objectiveFunctionGeometry(x):

    global objective_function_string

    global ground_truth, standardised_roi_ground_truth
    global best_fitness, best_fitness_id, fitness_function_call_id, ↴
    ↪evolution_fitness, evolution_parameters

    objectives = []

    for ind in x:
        xray_image = updateXRayImage(ind)
        corrected_xray_image = applyLogScaleAndNegative(xray_image[y_min_id: ↴
        ↪y_max_id, x_min_id:x_max_id])
        corrected_xray_image = corrected_xray_image.astype(np.single)

        if corrected_xray_image.min() != corrected_xray_image.max():
            standardised_corrected_xray_image = ↴
            ↪standardisation(corrected_xray_image)
        else:
            standardised_corrected_xray_image = corrected_xray_image

        # gX = cv2.Sobel(corrected_xray_image, ddepth=cv2.CV_32F, dx=1, dy=0)
        # gY = cv2.Sobel(corrected_xray_image, ddepth=cv2.CV_32F, dx=0, dy=1)
        # gX = cv2.convertScaleAbs(gX)
        # gY = cv2.convertScaleAbs(gY)
        # test_grad_magn = cv2.addWeighted(gX, 0.5, gY, 0.5, 0)

    ref_image = standardised_roi_ground_truth
    test_image = standardised_corrected_xray_image

    # ref_image = ref_grad_magn
    # if test_grad_magn.min() != test_grad_magn.max():
    #     test_grad_magn = standardisation(test_grad_magn)
    # else:
    #     test_image = test_grad_magn

    zncc = np.mean(standardised_roi_ground_truth * ↴
    ↪standardised_corrected_xray_image)
    dzncc = (1.0 - zncc) / 2.0
    objective = dzncc
```

```

#mae = np.mean(np.abs(standardised_roi_ground_truth - standardised_corrected_xray_image))
#objective = mae

#rmse = math.sqrt(np.mean(np.square(standardised_roi_ground_truth - standardised_corrected_xray_image)))
#objective = rmse

#ssim_value = ssim(standardised_roi_ground_truth, standardised_corrected_xray_image, data_range=standardised_roi_ground_truth.max() - standardised_roi_ground_truth.min())
#dssim = (1.0 - ssim_value) / 2.0
#objective = dssim

# # Avoid div by 0
# offset1 = min(standardised_roi_ground_truth.min(), standardised_corrected_xray_image.min())
# offset2 = 0.01 * (standardised_roi_ground_truth.max() - standardised_roi_ground_truth.min())
# offset = offset2 - offset1

#mape_value = mape(standardised_roi_ground_truth + offset, standardised_corrected_xray_image + offset)
#objective = mape_value

#mi = normalized_mutual_information(standardised_roi_ground_truth, standardised_corrected_xray_image)
#dmi = (1.0 - mi) / 2.0
#objective = dmi

objectives.append(objective)

return objectives

```

```

[56]: class SingleObjectiveProblem(Problem):

    def __init__(self):
        super().__init__(n_var=len(x_init), n_obj=1, n_constr=0, n_ieq_constr=0,
                        xl=xl, xu=xu, vtype=float)

    def _evaluate(self, x, out, *args, **kwargs):
        out["F"] = np.array(objectiveFunctionGeometry(x))
        out["F"].shape = [len(x), 1]

```

```
[57]: resetToDefaultParameters()
if os.path.exists(output_path + "/lungman-res-cmaes-X.dat") and os.path.
    exists(output_path + "/lungman-res-cmaes-F.dat"):
    res_cmaes_X = np.loadtxt(output_path + "/lungman-res-cmaes-X.dat")
    res_cmaes_F = np.loadtxt(output_path + "/lungman-res-cmaes-F.dat")
else:
    problem = SingleObjectiveProblem()


    r = np.array(xu) - np.array(xl);
    # #x0 = denormalize(np.random.random(problem.n_var), problem.xl, problem.xu)
    x0 = (np.random.random(problem.n_var) * r) + np.array(xl)
    print(x0)
    # #x0 = np.zeros(problem.n_var)
    algorithm = CMAES(x0=x0)

    res = minimize(problem,
                  algorithm,
                  ('n_iter', 1000),
                  seed=1,
                  save_history=True,
                  verbose=True)

    res_cmaes_X = res.X
    res_cmaes_F = res.F

    np.savetxt(output_path + "/lungman-res-cmaes-X.dat", res_cmaes_X)
    np.savetxt(output_path + "/lungman-res-cmaes-F.dat", res_cmaes_F)
```

```
[58]: xray_image_rmse_cmaes = applyLogScaleAndNegative(updateXRayImage(res_cmaes_X))
displayLinearPowerScales(xray_image_rmse_cmaes, "f", "f")
```



```
[59]: # Apply the transformation
xray_image_rmse_cmaes = applyLogScaleAndNegative(updateXRayImage(res_cmaes_X,
                                                               restore_transformation=False));

[60]: standardised_corrected_xray_image = standardisation(xray_image_rmse_cmaes[y_min_id:y_max_id, x_min_id:x_max_id])

imwrite(output_path + "/standardised_corrected_xray_image.tif", standardised_corrected_xray_image.astype(np.single))
hist_ref = np.histogram(standardised_roi_ground_truth, 100)[0]

standardised_roi_ground_truth = standardised_roi_ground_truth.astype(np.single)
gX = cv2.Sobel(standardised_roi_ground_truth, ddepth=cv2.CV_32F, dx=1, dy=0)
gY = cv2.Sobel(standardised_roi_ground_truth, ddepth=cv2.CV_32F, dx=0, dy=1)
gX = cv2.convertScaleAbs(gX)
gY = cv2.convertScaleAbs(gY)
ref_grad_magn = cv2.addWeighted(gX, 0.5, gY, 0.5, 0)

hist_ref = np.histogram(ref_grad_magn, 100)[0]

[61]: gvxr.displayScene()
gvxr.useNegative()

gvxr.setZoom(1339.6787109375)
gvxr.setSceneRotationMatrix([0.8227577805519104, 0.1368587613105774, -0.
                           →5516625642776489, 0.0, -0.5680444240570068, 0.23148967325687408, -0.
                           →7897683382034302, 0.0, 0.01961756870150566, 0.9631487131118774, 0.
                           →26820749044418335, 0.0, 0.0, 0.0, 0.0, 1.0])

gvxr.setWindowBackGroundColour(0.5, 0.5, 0.5)

gvxr.displayScene()

[62]: # gvxr.renderLoop()

[63]: # print(gvxr.getZoom())
# print(gvxr.getSceneRotationMatrix())

[64]: screenshot = (255 * np.array(gvxr.takeScreenshot())).astype(np.uint8)


[65]: plt.figure(figsize=(10,10))
plt.title("Screenshot")
plt.imshow(screenshot)
plt.axis('off')
```

```

plt.tight_layout()

plt.savefig(output_path + '/lungman-optimised-screenshot-beam-on.pdf')
plt.savefig(output_path + '/lungman-optimised-screenshot-beam-on.png')


```



```

[66]: fullCompareImages(standardised_roi_ground_truth,
                      standardised_corrected_xray_image,
                      "gVirtualXRay",
                      # output_path + "/lungman-projection-rmse-cmaes", ↴
                      ↴vmin=standardised_roi_ground_truth.min(), vmax=standardised_roi_ground_truth.
                      ↴max(),
                      output_path + "/lungman-projection-rmse-cmaes", ↴
                      ↴detector_element_spacing, vmin=-1.5, vmax=1.5)

```



```
[67]: ref_diag = np.diag(standardised_roi_ground_truth)
test_diag = np.diag(standardised_corrected_xray_image)

plt.figure(figsize=(15, 7))

ax = plt.subplot(111)


ax.set_title("Diagonal profiles")

ax.plot(ref_diag, label="Projection from ")
ax.plot(test_diag, label="gVirtualXRay (without post-processing)")

ax.legend(loc='best',
           ncol=1, fancybox=True, shadow=True)

# plt.legend()

plt.savefig(output_path + '/lungman-profiles-projection-rmse-cmaes.pdf')
plt.savefig(output_path + '/lungman-profiles-projection-rmse-cmaes.png')
```



8 Compute the magnification

$magnification = \frac{SID}{SOD}$ with SID the source to imager distance and SOD the source to object distance.

```
[68]: source_position = gvxr.getSourcePosition("mm")
detector_position = gvxr.getDetectorPosition("mm")

object_bbox = gvxr.getNodeAndChildrenBoundingBox("root", "mm")
object_position = [(object_bbox[0] + object_bbox[3]) / 2,
                   (object_bbox[1] + object_bbox[4]) / 2,
                   (object_bbox[2] + object_bbox[5]) / 2
                  ]

source_imager_distance = distance.euclidean(source_position, detector_position)
source_object_distance = distance.euclidean(source_position, object_position)

magnification = source_imager_distance / source_object_distance
```

```
[69]: print("SID:", source_imager_distance, "mm")
print("SOD:", source_object_distance, "mm")
print("magnification:", magnification)
```

SID: 1662.6100648426484 mm
 SOD: 1541.0320332447461 mm
 magnification: 1.0788939029008446

9 Compute the pixel size in the object plane

```
[70]: spacing = source_imager_distance * detector_element_spacing / ↴
       source_object_distance

print("spacing in the object plane (in mm):", spacing)
```

spacing in the object plane (in mm): [0.2157788 0.2157788]

10 Post-processing using image sharpening

We can see from the real image that an image sharpening filter was applied. We will implement one and optimise its parameters.

```
[71]: def sharpen(image, ksize, shift, scale):
        details = image - gaussian(image, ksize)

        return (details + shift) * scale
```

10.1 Define an objective function

```
[72]: # Compute the X-ray image
xray_image = np.flip(gvxx.computeXRayImage())

# Flat-field
xray_image_flat = xray_image[y_min_id:y_max_id, x_min_id:x_max_id] /_
    ↳total_energy_in_MeV
```

```
[73]: def objectiveFunctionSharpen(x):

    global xray_image_flat
    global standardised_roi_ground_truth

    objectives = []

    for parameters in x:

        # Retrieve the parameters
        shift1, scale1, sigma1, sigma2, shift2, scale2 = parameters

        # Process the image
        contrast = (xray_image_flat + shift1) * scale1
        details = sharpen(xray_image_flat, (sigma1, sigma2), shift2, scale2)
        sharpened = contrast + details
        sharpened[sharpened < 1e-9] = 1e-9
        log_image = np.log(contrast + details)
        negative = log_image * -1
        normalised = standardisation(negative)

        # Return the objective
        objective = math.sqrt(mean_squared_error(standardised_roi_ground_truth,_
            ↳normalised))
        # objective = math.sqrt(mean_squared_error(ref_grad_magn,_
        ↳test_grad_magn))
        # objective = math.sqrt(mean_squared_error(hist_ref, hist_test))

        # objective = ref_grad_magn.sum() - test_grad_magn.sum()
        # objective *= objective
        # objective *= -1

        # mi = normalized_mutual_info_score(hist_ref, hist_test)
        # dmi = (1.0 - mi) / 2.0
        # objective = dmi
        objectives.append(objective)

    return objectives
```

10.2 Minimise the objective function

```
[74]: sigma1 = 2
sigma2 = 2
alpha = 10.5
shift = 0
scale = 1

xl = [0, 0, 0, -5, 0]
xu = [10, 10, 15, 5, 2]
x_init = [sigma1, sigma2, alpha, shift, scale]

#0.003937458431052107 4012.600582499311 3.916470602237863 0.28374201341640476 6.
→851503972136956 6.8658686847669045e-09

xl = [0, 1e-9, 1, 1, 0, 0.5]
xu = [10000, 10000, 10, 10, 10000, 10000]
x_init = [0, 1, 3, 3, 0, 1]
```

```
[75]: class SingleObjectiveProblem(Problem):

    def __init__(self):
        super().__init__(n_var=len(x_init), n_obj=1, n_constr=0, n_ieq_constr=0, u
→xl=xl, xu=xu, vtype=float)
        # super().__init__(n_var=len(x_init), n_obj=1, n_constr=0, u
→n_ieq_constr=0, vtype=float)

    def _evaluate(self, x, out, *args, **kwargs):

        out["F"] = np.array(objectiveFunctionSharpen(x))
        out["F"].shape = [len(x), 1]
```

```
[76]: # The registration has already been performed. Load the results.
if os.path.isfile(output_path + "/lungman_postprocess.dat"):

    shift1, scale1, sigma1, sigma2, shift2, scale2 = np.loadtxt(output_path + "/"
→lungman_postprocess.dat")

else:
    # Optimise
    # timeout_in_sec = 20 * 60 # 20 minutes
    # opts = cma.CMAOptions()
    # opts.set('tolfun', 1e-10)
    # opts['tolx'] = 1e-10
    # opts['timeout'] = timeout_in_sec
    # opts['bounds'] = [xl, xu]
```

```

# opts['CMA_stds'] = []

# for min_val, max_val in zip(opts['bounds'][0], opts['bounds'][1]):
#     opts['CMA_stds'].append(abs(max_val - min_val) * 0.05)

problem = SingleObjectiveProblem()
x0 = denormalize(np.random.random(problem.n_var), problem.xl, problem.xu)

algorithm = CMAES(x0=x0,
                   sigma=0.5,
                   restarts=2,
                   maxfevals=np.inf,
                   tolfun=1e-6,
                   tolx=1e-6,
                   restart_from_best=True,
                   bipop=True)

res = minimize(problem,
               algorithm,
               ('n_iter', 100),
               seed=1,
               save_history=True,
               verbose=True)

# # Optimise
# res = minimize(problem,
#                 algorithm,
#                 ('n_iter', 100),
#                 # seed=1,
#                 save_history=False,
#                 verbose=True)

print(f"Best solution found: \nX = {res.X}\nF = {res.F}\nCV= {res.CV}")

# Save the parameters
shift1, scale1, sigma1, sigma2, shift2, scale2 = res.X
np.savetxt(output_path + "/lungman_postprocess.dat", [shift1, scale1, sigma1, sigma2, shift2, scale2], header='shift1, scale1, sigma1, sigma2, shift2, scale2')

```

10.3 Apply the result of the optimisation

```
[77]: print(shift1, scale1, sigma1, sigma2, shift2, scale2) #      return scale * (shift1  
    ↪+ image) + alpha * details
```

```
# Process the image  
contrast = (xray_image_flat + shift1) * scale1  
details = sharpen(xray_image_flat, (sigma1, sigma2), shift2, scale2)  
sharpened = contrast + details  
sharpened[sharpened < 1e-9] = 1e-9  
log_image = np.log(contrast + details)  
negative = log_image * -1  
normalised = standardisation(negative)
```

```
0.0038399819169387573 9232.063565472905 7.339600967177158 7.924408749898391  
11.289006641605752 0.6185266797462592
```

```
[78]: font = {'size' : 12.5}  
        }  
matplotlib.rcParams['font', **font]
```

```
[79]: new_image_without_post_process = standardised_corrected_xray_image *  
    ↪raw_reference.std()  
new_image_without_post_process -= standardised_corrected_xray_image.mean()  
new_image_without_post_process += raw_reference.mean()  
  
new_image_with_post_process = normalised * raw_reference.std()  
new_image_with_post_process -= normalised.mean()  
new_image_with_post_process += raw_reference.mean()
```

```
[80]: old_zncc = np.mean(standardised_roi_ground_truth *  
    ↪standardised_corrected_xray_image)  
new_zncc = np.mean(standardised_roi_ground_truth * normalised)
```

```
[81]: old_mape = mape(raw_reference, new_image_without_post_process)  
new_mape = mape(raw_reference, new_image_with_post_process)
```

```
[82]: old_ssim = ssim(raw_reference, new_image_without_post_process,  
    ↪data_range=raw_reference.max() - raw_reference.min())  
new_ssim = ssim(raw_reference, new_image_with_post_process,  
    ↪data_range=raw_reference.max() - raw_reference.min())
```

```
[83]: print("ZNCC before sharpening:", str(100 * old_zncc) + "%")  
print("ZNCC after sharpening:", str(100 * new_zncc) + "%")  
print()
```

```
print("MAPE before sharpening:", str(100 * old_mape) + "%")
```

```

print("MAPE after sharpening:", str(100 * new_mape) + "%")
print()

print("SSIM before sharpening:", str(old_ssim))
print("SSIM after sharpening:", str(new_ssim))

```

ZNCC before sharpening: 98.33418136180634%
ZNCC after sharpening: 98.88902710233528%

MAPE before sharpening: 1.7749609879510588%
MAPE after sharpening: 1.571659798523078%

SSIM before sharpening: 0.9394900050793231
SSIM after sharpening: 0.940002487387807

```

[84]: ground_truth_diag = np.diag(raw_reference)
gvxr_diag = np.diag(new_image_with_post_process)

print(ground_truth_diag.shape)
x = np.linspace(0, ground_truth_diag.shape[0], ground_truth_diag.shape[0]) * ↴
    spacing[0]
plt.figure(figsize=(15, 5))

ax = plt.subplot(111)

ax.set_title("Diagonal profiles")

ax.plot(x, ground_truth_diag, label="Ground truth (Digital Radiography)")
ax.plot(x, gvxr_diag, label="gVirtualXRay")


ax.legend(loc='best',
          ncol=1, fancybox=True, shadow=True)

plt.xlabel("Pixel position (in mm)")
plt.ylabel("Pixel intensity")

plt.savefig(output_path + '/lungman-profiles-projection-postprocessing.pdf')
plt.savefig(output_path + '/lungman-profiles-projection-postprocessing.png')

```

(1511,)



```
[85]: window_centre = int(volume.GetMetaData("0028|1050").split("\\\\")[1]) # Use 0 for normal, 1 for harder, 2 for softer
window_width = int(volume.GetMetaData("0028|1051").split("\\\\")[1]) # Use 0 for normal, 1 for harder, 2 for softer

print("Window Center used: ", window_centre)
print("Window Width used: ", window_width)


vmin = window_centre - window_width / 2
vmax = window_centre + window_width / 2

view_position = volume.GetMetaData("0018|5101")
```

Window Center used: 6032

Window Width used: 2245


```
[86]: fullCompareImages(raw_reference,
                      new_image_with_post_process,
                      "gVirtualXRay",
                      output_path + "/lungman-projection-harder", spacing,
                      log=False, vmin=vmin, vmax=vmax)
```



```
[87]: plt.figure(figsize= (20,10))
xrange=range(raw_reference.shape[1])
yrange=range(raw_reference.shape[0])

plt.xlabel("Pixel position\n(in mm)")
plt.ylabel("Pixel position\n(in mm)")
plt.title(view_position + " projection\nkVp:" + str(kVp))
plt.imshow(raw_reference, cmap="gray",
           vmin=vmin, vmax=vmax,
           extent=[0,(raw_reference.shape[1]-1)*spacing[0],0,(raw_reference.
           ↵shape[0]-1)*spacing[1]])
plt.colorbar(orientation='vertical')

plt.savefig(output_path + '/lungman_experimental_DX_image-harder.pdf')
plt.savefig(output_path + '/lungman_experimental_DX_image-harder.png')
```



```
[88]: plt.figure(figsize= (20,10))
xrange=range(raw_reference.shape[1])
```


```

yrange=range(raw_reference.shape[0])

plt.xlabel("Pixel position\n(in mm)")
plt.ylabel("Pixel position\n(in mm)")
plt.title(view_position + " projection\nkVp:" + str(kVp))
plt.imshow(new_image_with_post_process, cmap="gray",
           vmin=vmin, vmax=vmax,
           extent=[0,(raw_reference.shape[1]-1)*spacing[0],0,(raw_reference.
           ↵shape[0]-1)*spacing[1]])
plt.colorbar(orientation='vertical')

plt.savefig(output_path + '/lungman_simulated_DX_image-harder.pdf')
plt.savefig(output_path + '/lungman_simulated_DX_image-harder.png')

```



```

[89]: runtimes = []

resetToDefaultParameters()
setTransformations(res_cmaes_X)

```

```

for i in range(25):
    start_time = datetime.datetime.now()

    raw_x_ray_image = np.array(gvxr.computeXRayImage())

    # Process the image
    xray_image_flat = raw_x_ray_image / total_energy_in_MeV
    contrast = (xray_image_flat + shift1) * scale1
    details = sharpen(xray_image_flat, (sigma1, sigma2), shift2, scale2)
    sharpened = contrast + details
    sharpened[sharpened < 1e-9] = 1e-9
    log_image = np.log(contrast + details)
    negative = log_image * -1

    end_time = datetime.datetime.now()
    delta_time = end_time - start_time
    runtimes.append(delta_time.total_seconds() * 1000)

```

[90]:

```

runtime_avg = round(np.mean(runtimes))
runtime_std = round(np.std(runtimes))

raw_x_ray_image = np.array(raw_x_ray_image)

```

[91]:

```

print("Lungman PA view & Real digital radiograph & " +
      "{0:0.2f}".format(100 * new_mape) + "\\\%" & " " +
      "{0:0.2f}".format(100 * new_zncc) + "\\\%" & " " +
      "{0:0.2f}".format(new_ssim) + " " & "$" +
      str(raw_x_ray_image.shape[1]) + " \\\times " + str(raw_x_ray_image.
      ↪shape[0]) + "$ " & " N/A " &" +
      str(number_of_triangles) + " " & " " +
      "$" + str(runtime_avg) + " \\\pm " + str(runtime_std) + "$ & N/A \\\\"")

```

Lungman PA view & Real digital radiograph & 1.57\% & 98.89\% & 0.94
& \$1871 \times 1881\\$ & N/A & 20973540 & \$518 \pm 60\\$ & N/A \\

11 Visualise the virtual patient


[92]:

```

fname = output_path + "/visualisation.png"

if os.path.exists(fname):
    display(Image(fname, width=1000))
else:
    plot = visualise(use_log=True)
    plot.display()

```



```
[93]: if not os.path.exists(fname):
    if plot is not None:
        plot.fetch_screenshot()

    data = base64.b64decode(plot.screenshot)
    with open(fname, 'wb') as fp:
        fp.write(data)
```

11.1 All done

Destroy the window

```
[94]: gvxr.destroyAllWindows()
```